

SURCAM

Dental Pro.

CATALOG



SURCAM
Dental Pro.

CATALOG

Foreword	Dental Implants [Explanation]	Dental Implants [Products]	Prosthetic Introduction
06-17	18-23	24-39	36-59
<hr/> History 06 Background 08 Executives 10 Management 12 Production / Quality Management 13 Customer Service 14 Specifications 16 Prosthetic Products 17	<hr/> The S-Type 20 The C-Type 21 Surface Treatment Configuration 22 Implant Packaging 23	<hr/> The S-Type 28 Drilling Procedure [S-Type] 30 The C-Type 32 Drilling Procedure [C-Type] 34	<hr/> Drills + Sets 36 S.O.S. Set 37 SURCAM's Advantages 38 Closed Tray Transfer 39 Wise Click 40 Multi System Concept 48 Multi-Unit Restoration 49 Holder Bar 56


Internal Hex [Prosthetic Items]	Wrenches & Components	Pro.Positions	Internal Three Leg Fit [Prosthetic Items]	External Hex [Prosthetic Items]
60-91	92-101	102-113	114-135	136-147
Wise Click 64 Multi System Concept 66 Multi-Unit Restoration 68 Ball Attachment 72 Temporary Abutments 75 Titanium Abutments 76 Transfer 83 Casting Abutments 85 Analog 86 3D Scanning TiBase 87 Healing Cap 89 Prosthetic Screws 91	Manual Tools 96 Machine Tools 97 Ratchet Wrench 98 Crocks-Key 99 Adapter 100	Pro.Positions [Explanation] 104 Pro.Positions [Products] 107	Inverse Angular Abutments 116 Wise Click 120 Multi System Concept 122 Multi-Unit Restoration 126 Ball Attachment 128 Aesthetic Abutments 130 Casting Abutments 132 Analog 133 3D Scanning TiBase 134 Healing Cap 135	Wise Click 140 Multi-Unit Restoration 141 Ball Attachment 142 Abutments 143 Casting Abutments 144 Analog 145 Healing Cap 146 Transfer 147



SURCAM is a fast moving, dynamic company with a focus on designing, developing and manufacturing high quality medical instrumentation. We specialize in the manufacturing of medical equipment and titanium abutments for dental implants.

SURCAM was established in 2002 by Mr. Asher Ben Shabbat, who had a desire to create a place where medical experts and clinicians can realize their ideas and medical innovations. Mr. Ben Shabbat offered his own vast mechanical knowledge and many years of experience in the development and production of medical equipment and created SURCAM, which is in essence a high quality environment for the development of outstanding mechanical devices for medical applications.


SURCAM
DENTAL

BA  KGROUND

SURCAM's background is rooted in innovation and relies on great experience in the Medical Device world.

SURCAM takes pride in a rich tradition of innovation as it provides high-quality products with simplified procedures at value-added prices.

SURCAM DENTAL started to operate in the field of dental implant's industry in 2006 with the launch of its first prosthetic system. thanks to its sophisticated manufacturing processes, it continues to develop and thrive. The company has a state-of-the-art manufacturing site that is both technologically advanced and cost efficient. This facility, houses CNC machines that are only configured to produce families of parts (rather than individual parts), they meet demanding tolerances while operating in a lights-out manufacturing environment – thus achieving extraordinary efficiency and quality.





Asher Ben Shabat - FOUNDER and CEO - R&D, PRODUCTION

With over 30 years of technical R&D and business management experience in the Medical Devices Industry, and extensive and valuable experience in medical equipment development and production for major medical producers around the globe, Mr. Ben Shabbat offers his knowhow and understanding of medical procedures to surgeons and technical experts and helps them with the most difficult task: transforming ideas and sketches into prototypes and into the actual production of high quality medical devices.



Gil Lavi - PARTNER / V.P - MARKETING, LOGISTIC & QA

Offers a welth of 20 years of Marketing, sales and business management experience in a wide range of professional customer markets.



Management

SURCAM is managed by Mr. Ben Shabbat and an experienced staff. We are alert to the latest improvements in forming technologies and are equipped with the most recent, up to date modern equipment that is capable of manufacturing any desired part/shape, by utilizing the most efficient and cost effective methods.

SURCAM manufacturing techniques are constantly evolving to meet future challenges of the industry's next generation requirements at a competitive price.

SURCAM prides itself of producing a wide range of high quality Dental Abutments, while thoroughly understanding the actual dental procedures.

SURCAM offers the biggest range of Abutment designs that fit major dental implant manufactures, Dental Abutments, Screws, devices, tools and fixtures.

Production


All of SURCAM's Internal Hexed full line of Implants, Prosthetic Components and Abutments are compatible with standard platforms: 3.75, 4.2, 4.5mm [Zimmer, MIS, AlphaBio, BioHorizen].SURCAM's manufacturing techniques are constantly evolving to the future challenges of the industry's next generation requirements at a competitive price. All of our products are manufactured in-house.

Quality Management

SURCAM MEDICAL DEVICES AND DEVELOPMENT LTD Quality System is compliant and certified under the following international Standards:

ISO 13485:2003
CE 0473





CUSTOMER SERVICE

One of SURCAM's most important founding principles is customer service and support.

This means fulfilling individual customer needs. This is accomplished by establishing flexible operations systems which allows products to be easily tailored to specific needs while still adhering to short delivery times. It also means providing a quick response to inquiries - whether technical or commercial in nature. We pride ourselves of precise accurate answers within minutes or hours.

We are constantly working to improve our already high standards of quality, performance and reliability, which are backed by superb service and support.

SPECIFICATIONS


- Dental Implants titanium grade 5 with SLA surface treatment.
- Aesthetic Abutment body functions as impression coping.
- All SURCAM products Includes the Prosthetic Screw.
- Installation using large external hex drivers (1.3mm).
- Non-rotational design for standard, internal-hex implants.
- Compatible with internal hex platform: 3.75, 4.2 & 4.5 mm Implants.
[Zimmer®, MIS®, AlphaBio®, A.B Dental®, Bio Horizon®].
- Compatible - Internal Three Leg Fit platform: 3.5, 4.3 & 5.0 mm Implants.
[Noble Biocare®, Replace Select®, Tri-lob®].
- Compatible with external hex platform: 4.1 & 5.0 mm Implants.
[Brenemark®].
- Compatible with internal hex conical connection. platform: 3.5 & 4.0 mm Implants.
[Ti Design®].
- Made of Commercial grade 5 Titanium alloy according to ASTM-F-136 02A COND A.
- Stainless Steel AISI 303 alloy according to ASTM A582-95B COND A.

Advantages and Priorities

In order to preserve the original structure of the titanium alloy, SURCAM utilizes advanced CAD/CAM methods in the production process.

These methods effectively reduce undesirable deformations of the alloy structure.


PROSTHETIC PRODUCTS



All of SURCAM's full line of internal connections prosthetic components and abutment designs are compatible with: **Internal Hex -3.75 & 4.5 platform, Internal Three Leg Fit -3.5 , 4.3 & 5.0 platform, and External Hex - 4.1 & 5.0 platform.**

SURCAM offers the restorative dentist the convenience of using a wide range of prosthetic platforms when various types of abutments are required. Each abutment design is ideally engineered to suit a variety of surgical cases.


SURCAM's philosophy is that properly matching the abutment to its final shape in the machining stage is one of the most important elements of the entire restoration procedure.



Dental
MPLANTS

The S-Type Spiral body is designed with self-drilling capability for high primary stability.

The SURCAM self-tapping S-Type implants orientation can be adjusted during insertion with high initial stability even in soft bone, which enables a smaller osteotomy, resulting in minimal bone loss and reduced trauma. The S-type implants are the optimal solution for immediate loading restoration.



	A	B
S33008INX	3.30	2.92
S33010INX	3.30	2.92
S330115INX	3.30	2.92
S33013INX	3.30	2.92

	A	B
S37508INX	3.78	2.01
S37510INX	3.78	2.01
S375115INX	3.78	2.01
S37513INX	3.78	2.01
S37516INX	3.78	2.01
S37518INX	3.78	2.01


	A	B
S42008INX	4.17	2.40
S42010INX	4.17	2.40
S420115INX	4.17	2.40
S42013INX	4.17	2.40
S42016INX	4.17	2.40
S42018INX	4.17	2.40

	A	B
S50008INX	4.98	3.00
S50010INX	4.98	3.00
S500115INX	4.98	3.00
S50013INX	4.98	3.00
S50016INX	4.98	3.00
S50018INX	4.98	3.00

The C-Type

Classic Thread Parallel shape.

The C-Type implant's cylindrical body is designed to maximize initial stability in soft bone. The shallow cutting edges on its threads are specifically designed for soft bone indications and minimize pressure on hard bone.



	A	B
C33008INX	3.30	2.05
C33010INX	3.30	2.05
C330115INX	3.30	2.05
C33013INX	3.30	2.05
C33016INX	3.30	2.05


	A	B
C37508INX	3.78	2.65
C37510INX	3.78	2.65
C375115INX	3.78	2.65
C37513INX	3.78	2.65
C37516INX	3.78	2.65

	A	B
C42008INX	4.19	3.13
C42010INX	4.19	3.13
C420115INX	4.19	3.13
C42013INX	4.19	3.13
C42016INX	4.19	3.13


	A	B
C50008INX	4.98	3.60
C50010INX	4.98	3.60
C500115INX	4.98	3.60
C50013INX	4.98	3.60
C50016INX	4.98	3.60

SURFACE TREATMENT CONFIGURATION

SURCAM's surface treatment of Dental Implants is at a very high level, and includes processes such as dry alumina blasting and double acid etching (SLA type of dental implants surface treatment).




The surface of the implant at low magnification (X50) in back scattering electrons.










The Implant's surface at high magnification (X2000) in secondary electrons.

IMPLANT PACKAGING




1 Single Implant box - Implant type Label [S-Type / U-Type].

2 10 Package box - your Implant type is marked.




A teal-colored geometric shape, resembling a parallelogram or a trapezoid, is positioned behind the text "DENTAL".

DENTAL IMPLANTS



Our Implants Variety









S-TYPE

Spiral body design with self-drilling capability for high primary stability

S-TYPE CONICAL IMPLANT

Dia. 3.30 mm ●

			
S33008INX Dia. = 3.30 mm L = 8.0 mm	S33010INX Dia. = 3.30 mm L = 10.0 mm	S330115INX Dia. = 3.30 mm L = 11.5 mm	S33013INX Dia. = 3.30 mm L = 13.0 mm

S-TYPE CONICAL IMPLANT

Dia. 3.75 mm ●

					
S37508INX Dia. = 3.75 mm L = 8.0 mm	S37510INX Dia. = 3.75 mm L = 10.0 mm	S375115INX Dia. = 3.75 mm L = 11.5 mm	S37513INX Dia. = 3.75 mm L = 13.0 mm	S37516INX Dia. = 3.75 mm L = 16.0 mm	S37518INX Dia. = 3.75 mm L = 18.0 mm

S-TYPE CONICAL IMPLANT

Dia. 4.20 mm ●

					
S42008INX Dia. = 4.2 mm L = 8.0 mm	S42010INX Dia. = 4.2 mm L = 10.0 mm	S420115INX Dia. = 4.2 mm L = 11.5 mm	S42013INX Dia. = 4.2 mm L = 13.0 mm	S42016INX Dia. = 4.2 mm L = 16.0 mm	S42018INX Dia. = 4.2 mm L = 18.0 mm

S-TYPE CONICAL IMPLANT

Dia. 5.0 mm ●


					
S50008INX Dia. = 5.0 mm L = 8.0 mm	S50010INX Dia. = 5.0 mm L = 10.0 mm	S500115INX Dia. = 5.0 mm L = 11.5 mm	S50013INX Dia. = 5.0 mm L = 13.0 mm	S50016INX Dia. = 5.0 mm L = 16.0 mm	S50018INX Dia. = 5.0 mm L = 18.0 mm



S-TYPE CONICAL IMPLANT

Dia. 6.0 mm ●

				
S60008INX Dia. = 6.0 mm L = 8.0 mm	S60010INX Dia. = 6.0 mm L = 10.0 mm	S600115INX Dia. = 6.0 mm L = 11.5 mm	S60013INX Dia. = 6.0 mm L = 13.0 mm	S60016INX Dia. = 6.0 mm L = 16.0 mm

Drilling Sequence - S-Type






- Drill to full depth of the Implant (*See illustrations 1 & 2)
- Drill to part depth of the Implant (*See illustrations 1 & 2)

▶ Drill only through the cortical bone, the drill should not be used to full depth.


The chart above illustrates the recommended drilling sequence for SURCAM's implants. Clinical judgement with respect to each individual patient must supersede any recommendation made.

illustrations 1



For bone type I II

illustrations 2



For bone type III IV

2.50 mm Depth
Small Dia. Drill

C-TYPE

Classic thread Parallel Type

C-TYPE **NARROW** PARALLEL IMPLANT

Dia. 3.30 mm ●

			
C33008INX Dia. = 3.3 mm L = 8.0 mm	C33010INX Dia. = 3.3 mm L = 10.0 mm	C330115INX Dia. = 3.3 mm L = 11.5 mm	C33013INX Dia. = 3.3 mm L = 13.0 mm

C-TYPE PARALLEL IMPLANT

Dia. 3.75 mm ●

				
C37508INX Dia. = 3.75 mm L = 8.0 mm	C37510INX Dia. = 3.75 mm L = 10.0 mm	C375115INX Dia. = 3.75 mm L = 11.5 mm	C37513INX Dia. = 3.75 mm L = 13.0 mm	C37516INX Dia. = 5.0 mm L = 16.0 mm


C-TYPE PARALLEL IMPLANT

Dia. 4.20 mm ●

				
C42008INX Dia. = 4.2 mm L = 8.0 mm	C42010INX Dia. = 4.2 mm L = 10.0 mm	C42115INX Dia. = 4.2 mm L = 11.5 mm	C42013INX Dia. = 4.2 mm L = 13.0 mm	C42016INX Dia. = 4.2 mm L = 16.0 mm

C-TYPE PARALLEL IMPLANT

Dia. 5.0 mm ●


				
C50008INX Dia. = 5.0 mm L = 8.0 mm	C50010INX Dia. = 5.0 mm L = 10.0 mm	C500115INX Dia. = 5.0 mm L = 11.5 mm	C50013INX Dia. = 5.0 mm L = 13.0 mm	C50016INX Dia. = 5.0 mm L = 16.0 mm


C-TYPE PARALLEL IMPLANT

Dia. 6.0 mm ●

		
C60008INX Dia. = 6.0 mm L = 8.0 mm	C60010INX Dia. = 6.0 mm L = 10.0 mm	C600115INX Dia. = 6.0 mm L = 11.5 mm

Drilling Sequence - C-Type







- Drill to full depth of the Implant (*See illustrations 1 & 2)
- Drill to part depth of the Implant (*See illustrations 1 & 2)

▶ Drill only through the cortical bone, the drill should not be used to full depth.

The chart above illustrates the recommended drilling sequence for SURCAM's implants. Clinical judgement with respect to each individual patient must supersede any recommendation made.



For **bone type II**



For **bone type III IV**

2.50 mm
Small Dia. Drill

Drills

Lance Drill

Dia. 1.50 Ø




DRL01515

Dia. = 1.50 Ø
L = 15.0 mm

Cylindrical Drill

Dia. 2.00 Ø



DR020130

Dia. = 2.00 Ø
L = 13.0 mm

Cylindrical Drill

Dia. 2.50 Ø



DR025130

Dia. = 2.50 Ø
L = 13.0 mm

Cylindrical Drill

Dia. 2.80 Ø



DR028130

Dia. = 2.80 Ø
L = 13.0 mm

Cylindrical Drill

Dia. 3.20 Ø



DR032130

Dia. = 3.20 Ø
L = 13.0 mm

Cylindrical Drill

Dia. 3.65 Ø



DR0365130

Dia. = 3.65 Ø
L = 13.0 mm

Cylindrical Drill

Dia. 4.10 Ø



DR041130

Dia. = 4.10 Ø
L = 13.0 mm

Cylindrical Drill

Dia. 4.30 Ø



DR043130

Dia. = 4.30 Ø
L = 13.0 mm

Cylindrical Drill

Dia. 5.20 Ø



DR052130


Dia. = 5.20 Ø
L = 13.0 mm

SURGICAL SET




SGK- 400106

KSH-1513 / KSH-2013 / MKM-1513 / MKS-0813 / KSH-1224 / KSH-1824
MIC-1224 / MIC-1824 / TRTS-02 / DR020115 / DR025115 / DR028115 / DR032115 DR0365115
/ TDE-400006 / DRM15016 / DR045115 DR041115 / DG-400103 / DGS-400104 / DGL-400105




Set


- 1 Implant with broken screw below the occlusal surface of the implants.




- 2 Place the drill guide (SG-172) in the implant hex.




- 3 Place the SOS Drill (SD-172) in the guide. Counter clockwise direction. Drilling speed: 700-1200Rpm. While removing the screw, use repeated tapping movements through the drill guide.



- 4 Place the SOS Climb Drill (SCD-172) in the guide. The speed must be under 80Rpm. Counter clockwise direction. While removing the screw, use repeated tapping movements through the drill guide.



- 5 Take the broken screw out!



S.O.S. Set - Screw Out Set


SURCAM's ADVANT/GES

1 Abutment central positioning.

3 The screw head is precisely designed for merging with the abutment shape.

2 Perfect match with all major implant companies Platforms.

4 The interface angle creates a bigger contact surface for perfect centering and repeatability when tightening the abutment screw in the lab or on the implant. The angle at the top of the screw creates a larger contact area which requires a high opening torque.



Closed Tray TRANSFER with Plastic Cap

1

External hex with conical lead

The hex is parallel to the implant hex, which enables better positioning achieving good aesthetic restoration.



3

Titanium body

The Titanium body can be used as a titanium abutment for temporary or final restoration.

4

5 Different heights



*

2

TiN coating

TiN coatings only modify surface properties and have no effect on the substrate properties or biomechanical functionality.

Biocompatibility:

Biocompatibility of TiN has been demonstrated in numerous published studies. It is known to improve friction, and has a significantly higher degree of surface scratch resistance compared to other materials.

5


Plastic Cap in two heights


The Plastic caps are machined in order to maintain the flexibility of the material and enable the release of the cap at up to 25° angle. It also maintains the accuracy.



6

Snap-on to secure positioning.






Overdenture Attachment

Instructions for use

- 1 The SURCAM system relies on a unique design of a hug ring made of PEEK, i.e. Poly Ether Ether Ketone (PEEK).
For more information go to:
<http://www.emdt.co.uk/article/use-peek-advanced-active-implants>.
- 2 The PEEK polymer shows excellent resistance to stress and Plastic deformations, which extends the life of the product beyond the lives of other alternatives products that are available on the market.
- 3 The unique mechanical and elastic properties of PEEK allow for a large diameter of spherical hug, resulting a very high grip force. The geometry of the titanium abutment (base unit ring) and the accompanying PEEK attachment create a large contact area, ensuring excellent stability of the prosthesis and comfort to the patients.



4

There is a variety of the attachments that can be used on the denture.
All the attachments connect to the Wise Click abutments:

**WCP-2348**

4.8 mm dia. PEEK attachment, to be mounted without metal housing in very narrow cases.

**WCH-4823**

4.8 mm dia. semi-protection of the PEEK attachment, secondary, preferred in narrow cases between the implants.

**WCH-5628**

5.6 mm dia. full-protection of the PEEK attachment. First choice to be preferred.

**WCP-2348-H**

HARD - Use upon request for better stability of the prosthesis.

**WCPF-2348-H**

POSITIONING - Use upon request for better stability of the prosthesis.


**WCPS-2348-H**

SOFT - For primary use.


- ☐ WISE CLICK STRAIGHT ABUTMENTS (Abutment bases connect to implants).



- ☐ WISE CLICK ANGULAR ABUTMENT (can be mounted at three different angles: 9°, 18°, 30°)





- ☐ WISE CLICK COVER (can be mounted at four different heights: 2, 3, 4, 5 mm)



- ☐ WISE CLICK ANGULAR ABUTMENT (can be mounted at three different angles: 9°, 18°, 30°)




■ ANGLE AND HEIGHT CLIP GAUGE (Use the Gauge Set for angle and height selection).




Related to pages: 71 + 76-77 + 80-81. [71 - Wise click, 76-77 - Multi-Unit, 80-81 - Ball attachment].


- 5 The wide contact area is possible because of the wide range of angles: 0°, 9°, 18°, 30° and wide range of heights: 1,2,3,4,5 & 6 mm. The availability of these options allows a near-perfect match for the contact areas between the titanium abutment (base unit ring) and the PEEK attachment.




- 6 When mounting the WISE CLICK SYSTEM, it is critical to verify up to 5° parallelism between the PEEK attachment and the base abutment.



- For placement of the Wise Click straight abutment use a standard 1.27 mm hex drive tool.




- **MULTI SYSTEM BASE CONNECTION** (Use the Plastic holder to place the base connection into the implant).



9° - 18° - 30°

See page # 52



1.0 mm - 4.0 mm

See page # 52

- **CLOSED TRAY IMPRESSION COPING**
Connect the transfer onto the abutment. Inject impression material and position the impression tray.



TWC-1105




AWC-1301

Send the impression tray to the laboratory with analogs.




WISE CLICK FIXATION

Protecting Nylon Sealed Ring (CAT NO: SCH-8540). should be mounted on the WISE CLICK housing slot, along with the WCPF-2348-H. **Use WCPF-2348-H only for positioning, not for attachment.**



Hold without pressure until material is set.





CONNECTION OF THE FINAL PROSTHETIC RESTORATION: Wise Click denture caps are placed. The Nylon Ring (SCH-8540) should be removed.

Note

Please be aware that the PEEK attachment should be replaced if exposed to a temperature higher than 100°C.

■ MAINTENANCE

- A Replacing the PEEKs holding attachment is very simple. Threading the WCT 900 special tool into the threaded hole in the PEEK attachment - will cause the ejection of the PEEK part, as shown.
- B Put the new peek attachment into the opposite side of the special tool.












































Note

In cases where there are significant differences in heights and angles, the patient should be instructed to press the prosthesis first at the front of mouth until hearing the engagement click and then tighten the back of the prosthesis until sensing the second engagement click.

MULTI SYSTEM CONCEPT

This system provides the clinician with large Adapters that are economical; A single Adapter can be used with a combination of up to 3 different families of over-denture sets.


 EMU-0910 9°			 EMU-1810 18°			 EMU-3010 30°		
Multi-Unit			Ball Attachment			Wise Click		
 MUC-0021	 MUC-0031		 BMA-0020	 BMA-0030		 WCC-0020	 WCC-0030	
 MUC-0041	 MUC-0051		 BMA-0040	 BMA-0050		 WCC-0040	 WCC-0050	
 PSU-0050	 PSU-0032	 TSU-1250	 BHM-5032	 BIS-0025		 WCH-4823	 WCP-2348-H	
 TMU-1000	 TMU-8005	 TMU-1200	 BHM-4730	 AMB-1201		 WCPF-2348-H	 WCPS-2348-H	
 AMU-1201	 TMU-0532	 HCC-4550				 SCH-8540	 WCP-2348	
 PSU-100015	 PSU-1050	 TSM-100013				 AWC-1301	 TWC-1105	
 CMU-100006	 CMU-100007							

Note Full description and information about the accessories find in pages # 72+77.


MULTI-UNIT RESTORATION

Screw Retained Restoration


- 1 Straight Multi-Unit - an abutment connection to an implant can be mounted at six different heights:




- 2 18°, 30° and 45° Angled Multi-Unit abutments are available in four different heights: H01, H02, H03, H04.






- 3 An Angular Adapter can be mounted at three different angles:




- 4 Multi-Unit Cover - Can be mounted at four different heights:



■ Accessories




- 5 Angle and height clip gauge - use the Gauge Set for angle and height selection.




Related to pages: 71 + 76-77 + 80-81. [71 - Wise click, 76-77 - Multi-Unit, 80-81 - Ball attachment].


- Using the angle gauges gives you an option to pre-plan the use of the abutment in order to achieve maximum parallelism.




- 6 Multi-Unit connection - use the holder to place the base connection into the implant.




- 7 Multi-Unit base connection - use the holder to place the base connection into the implant.



8 Match the Multi-Unit cover abutment height 2-5 mm as required.




9 Closed tray impression coping - with Plastic cap.




☐ Snap-on to secure positioning

☐ Use Closed-tray impression coping



- 10 Open tray impression coping - connect the transfer onto the abutment. Inject impression material and position the impression tray.




- Inject impression material and position the impression tray.
Send the impression tray to the laboratory.



- 11 In cases when a two-stage surgery is indicated:
Healing procedure - connect the temporary restoration healing cap.



- The original denture can be used on top of the healing caps while the temporary denture is being produced.



- 12 Connection of the temporary prosthetic restoration - verify tightening of 15 Ncm.




SMU-4023

Screw for
Multi-Unit
Sleeve

Note: The screw geometry was carefully designed and optimized to allow safe tightening and long service. Special attention must be paid while tightening the screw. Make sure to use an original SURCAM hex key with edges that do not show wear. It is recommended to purchase additional screws.


Simple
immediate
restoration
implantation



TSU-1250

Titanium
Sleeve


- 13 The sleeves height needs to be trimmed to the desired height.



Immediate Restoration Modular Passive HOLDER BAR

		
ISU-1201 H01	ISU-1202 H02	ISU-1203 H03


Sleeve - Titanium sleeve with three threads on the side walls at three different hights: H01, H02, H03



Cross section of the sleeve.
100°, 120°, 140° deviation.

		
ISU-1801 H01	ISU-1802 H02	ISU-1803 H03

Sleeve - Titanium sleeve with two threads on the side walls at three different hights: H01, H02, H03



Cross section of the sleeve.
180° deviation.



IBB-2530




IBM-3230




BBT-110

- > Ball Attachment to be mounted on the side wall of the sleeve.
- > Connection Bar to be modified to the space between the sleeves.
- > Bending Tool - For bending the Holder Bar [IBU-3330].

- 1 Place the Multi-unit on the Implants



- 2 Adjust the positions of the sleeves. Adjust sleeves mounted against each other




- 3 Attach Ball attachment to the sleeves mounted against each other



The titanium sleeve with Ball attachment mounted on the side walls using the standard 1.27 Hex tool.


4

The sleeve should be cut to the desired length. Modify the connecting bar so it fits into the space between the sleeves.




5

If necessary - the HOLDER BAR can be bent up to 5° in cold condition OR any angle needed up to 25°.




6

Attach the modified bar to the balls until you hear a "click" sound.



- 7 Fabricate the bridge on the modular passive frame for immediate restoration



- 8 Adjust the sleeves to the required height



G PRO GRIP



Internal hex

Compatible with internal hex platform: 3.75, 4.2 & 4.5 mm Implants.
[Zimmer®, MIS®, AlphaBio®, A.B Dental®].



Wise Click



Multi System
Concept



Multi-Unit
Restoration



Ball
Attachment

Components
for
Overdenture
Restoration



Titanium
Abutments



Temporary
Abutments

Components
for
Cement
Restoration



Casting
Abutments



Analog



3D Scanning
Interface

Components
for
Lab
Use



Transfers




Prosthetic
Screws



Healing
Cap

Complementary
Items



STRAIGHT WISE CLICK

Overdenture Attachment connection to implants

WISE CLICK MINOR - RETAINED OVERDENTURE ABUTMENTS

Abutment bases connection to implants

3.75 Platform

For Full Set add SPE to catalog number

					
WCM-0000 H = 0.3 mm	WCM-0010 H = 1.0 mm	WCM-0020 H = 2.0 mm	WCM-0030 H = 3.0 mm	WCM-0040 H = 4.0 mm	WCM-0050 H = 5.0 mm

WISE CLICK MINOR - SET PEEK

Abutment bases with peek attachment and housing covers

3.75 Platform

FULL SET

					
WCM-0000-SPE H = 0.3 mm	WCM-0010-SPE H = 1.0 mm	WCM-0020-SPE H = 2.0 mm	WCM-0030-SPE H = 3.0 mm	WCM-0040-SPE H = 4.0 mm	WCM-0050-SPE H = 5.0 mm

WISE CLICK - ACCESSORIES

Included in SET

						
WCH-5628 Stainless Steel Housing Cover	WCH-4823 Stainless Steel Housing	WCP-2348-H Hard Housing PEEK Attachment	WCPF-2348-H Positioning PEEK Attachment	WCPS-2348-H Soft Housing PEEK Attachment	WCP-2348 Direct PEEK Attachment	SCH-8540 Nylon Sealed Ring

Not Included in SET

		
TWC-1105 Transfer For Wise Click	WCT-900 Wise Click Tool	AWC-1301 Analog For Wise Click

ANGULAR WISE CLICK CONCEPT 9°, 18°, 30°

Wise Click Overdenture Attachment - SETS





9°

WISE CLICK ANGULAR - MULTI SYSTEM CONCEPT
Wise Click Angular overdenture attachment - Set

			
WCA-0920-M H02 = 3.2 mm 9° Wise Click Angular Multi System	WCA-0930-M H03 = 4.2 mm 9° Wise Click Angular Multi System	WCA-0940-M H04 = 5.2 mm 9° Wise Click Angular Multi System	WCA-0950-M H05 = 6.2 mm 9° Wise Click Angular Multi System




18°

WISE CLICK ANGULAR - MULTI SYSTEM CONCEPT
Wise Click Angular overdenture attachment - Set

			
WCA-1820-M H02 = 3.5 mm 18° Wise Click Angular Multi System	WCA-1830-M H03 = 4.5 mm 18° Wise Click Angular Multi System	WCA-1840-M H04 = 5.5 mm 18° Wise Click Angular Multi System	WCA-1850-M H05 = 6.5 mm 18° Wise Click Angular Multi System

30°

WISE CLICK ANGULAR - MULTI SYSTEM CONCEPT
Wise Click Angular overdenture attachment - Set

			
WCA-3020-M H02 = 4.5 mm 30° Wise Click Angular Multi System	WCA-3030-M H03 = 5.5 mm 30° Wise Click Angular Multi System	WCA-3040-M H04 = 6.2 mm 30° Wise Click Angular Multi System	WCA-3050-M H05 = 7.2 mm 30° Wise Click Angular Multi System

Designated information for Angular Multi System Adapters:

Internal Hex

3.75 Platform = WCA-....-M

4.5 Platform = WCA-....-W




ANGULAR MULTI SYSTEM CONCEPT 9°, 18°, 30°

Multi-Unit items Overdenture Attachments

AESTHETIC MINOR UNIT - SYSTEM ADAPTERS

Angular Multi System Adapter 9°, 18°, 30°




3.75 Platform

		
EMU-0910 9° Angle Multi System Adapter	EMU-1810 18° Angle Multi System Adapter	EMU-3010 30° Angle Multi System Adapter

AESTHETIC WIDE UNIT - SYSTEM ADAPTERS

Angular Multi System Adapter 9°, 18°, 30°


4.5 Platform

		
EWU-0910 9° Angle Multi System Adapter	EWU-1810 18° Angle Multi System Adapter	EWU-3010 30° Angle Multi System Adapter

GAUGE ANGLE AND HIGH MINOR

Clip Attachment

			
GAHM-00 Straight Gauge	GAHM-09 9° Angle Gauge	GAHM-18 18° Angle Gauge	GAHM-30 30° Angle Gauge


GAHM-SET Set includes the full range of Gauges





WISE CLICK CONNECTION

Wise Click connection for multi system concept

			
WCC-0020 2 mm Connection Multi System Wise Click	WCC-0030 3 mm Connection Multi System Wise Click	WCC-0040 4 mm Connection Multi System Wise Click	WCC-0050 5 mm Connection Multi System Wise Click





BALL MULTI ANGULAR

Ball Attachment connection for multi system concept

			
BMA-0020 2 mm Connection Multi System Ball	BMA-0030 3 mm Connection Multi System Ball	BMA-0040 4 mm Connection Multi System Ball	BMA-0050 5 mm Connection Multi System Ball

MULTI-UNIT CONNECTION

Screw-in connection for multi system concept

			
MUC-0021 2 mm Connection Multi System	MUC-0031 3 mm Connection Multi System	MUC-0041 4 mm Connection Multi System	MUC-0051 5 mm Connection Multi System


SMU-4023 Screw for Multi-Unit Sleeve

OVERDENTURE ABUTMENTS

For single tooth restoration

AESTHETIC CONNECTION ABUTMENTS

Screw-in abutments for single tooth restoration with hex locking for retention and positioning

		
ECA-6005 H = 0.5 mm	ECA-6015 H = 1.5 mm	ECA-6025 H = 2.5 mm

AESTHETIC PLASTIC SLEEVE

	
PES-1240 Plastic for WCA hex base	PES-1240-R Plastic for WCA round base

Related Screws:

SCS-3011	SCS-3012	SCS-3013
		

MULTI-UNIT STRAIGHT AESTHETIC ABUTMENTS

Abutments with smooth, non-locking surface for Multi-Unit Restoration

AESTHETIC MINOR UNIT

Straight Multi-Unit

3.75 Platform

					
EMU-0011 H = 1.0 mm	EMU-0021 H = 2.0 mm	EMU-0031 H = 3.0 mm	EMU-0041 H = 4.0 mm	EMU-0051 H = 5.0 mm	EMU-0061 H = 6.0 mm

AESTHETIC WIDE UNIT

Straight Multi-Unit

4.5 Platform

					
EWU-0011 H = 1.0 mm	EWU-0021 H = 2.0 mm	EWU-0031 H = 3.0 mm	EWU-0041 H = 4.0 mm	EWU-0051 H = 5.0 mm	EWU-0061 H = 6.0 mm


SMU-4023 Screw for Multi-Unit Sleeve

See the instructions Manual

ANGULAR MULTI SYSTEM CONCEPT 9°, 18°, 30°

Multi-Unit Overdenture attachment

9°

ANGULAR MULTI-UNIT SYSTEM CONCEPT

Aesthetic Multi-Unit Angular Overdenture Attachment - SET





			
EMA-0921-M H02 = 2.8 mm 9° Angular Multi System	EMA-0931-M H03 = 3.8 mm 9° Angular Multi System	EMA-0941-M H04 = 4.8 mm 9° Angular Multi System	EMA-0951-M H05 = 5.8 mm 9° Angular Multi System

18°

ANGULAR MULTI-UNIT SYSTEM CONCEPT

Aesthetic Multi-Unit Angular Overdenture Attachment

Mono-block

			
AMU-1801 H01 = 2.0 mm 18° Angular Multi-Unit	AMU-1811 H02 = 2.8 mm 18° Angular Multi-Unit	AMU-1821 H03 = 3.9 mm 18° Angular Multi-Unit	AMU-1831 H04 = 5.1 mm 18° Angular Multi-Unit

Two-Piece





			
EMA-1831-M H03 = 3.8 mm 18° Angular Multi System	EMA-1841-M H04 = 4.8 mm 18° Angular Multi System	EMA-1851-M H05 = 5.8 mm 18° Angular Multi System	EMA-1861-M H06 = 6.8 mm 18° Angular Multi System

30°

ANGULAR MULTI-UNIT SYSTEM CONCEPT

Aesthetic Multi-Unit Angular Overdenture Attachment

Mono-block

			
AMU-3001 H01 = 3.0 mm 30° Angular Multi-Unit	AMU-3011 H02 = 3.9 mm 30° Angular Multi-Unit	AMU-3021 H03 = 4.9 mm 30° Angular Multi-Unit	AMU-3031 H04 = 5.9 mm 30° Angular Multi-Unit

Two-Piece

			
EMA-3031-M H03 = 4.5 mm 30° Angular Multi System	EMA-3041-M H04 = 5.5 mm 30° Angular Multi System	EMA-3051-M H05 = 6.2 mm 30° Angular Multi System	EMA-3061-M H06 = 7.2 mm 30° Angular Multi System


Designated information for
Angular Multi System Adapters:**Internal Hex**

3.75 Platform = EMA-....-M

4.5 Platform = EMA-....-W

45°

ANGULAR MULTI-UNITAesthetic Multi-Unit Angular
Overdenture Attachment


AMU-100005 H01 = 4.5 mm 45° Angular Multi-Unit

MULTI-UNIT - ACCESSORIES

							
ISU-1201 H01	ISU-1202 H02	ISU-1203 H03	ISU-1801 H01	ISU-1802 H02	ISU-1803 H03	IBB-2530 Ball attachment	IBM-3230 Connection bar to be modified to the space between the sleeves
Ti Sleeve with 3 positioning options			Ti Sleeve with 2 positioning options				

							
PSU-0032 Plastic Sleeve screw included	TSU-1250 Titanium Sleeve screw included L = 12.0 mm	CMU-100006 Short slim titanium sleeve	CMU-100007 Long slim titanium sleeve L = 8.0 mm	PSU-1050 3D Scanning Lab device	PSU-100015 3D Scanning Intraoral device	PSU-0050 Parallel connection	TSM-100013 Titanium sleeve for welding

						
TMU-0532 Titanium Short Sleeve screw included	TMU-1000 Open Tray Short Transfer	TMU-1200 Open Tray Transfer screw included	TMU-8005 Closed-Tray Transfer	WQP-011 Plastic for Closed-Tray Transfer	HCC-4550 Healing Cap for Multi-Unit	AMU-1201 Analog for Multi-Unit

OVERDENTURE BALL ATTACHMENT

Connection to Implants

BALL MINOR STRAIGHT

3.75 Platform / Dia. 4.0 mm

					
BMS-4010 H = 1.0 mm Without Cap	BMS-4020 H = 2.0 mm Without Cap	BMS-4030 H = 3.0 mm Without Cap	BMS-4040 H = 4.0 mm Without Cap	BMS-4050 H = 5.0 mm Without Cap	BMS-4060 H = 6.0 mm Without Cap

BALL MINOR SET - NYLON

3.75 Platform / Dia. 4.0 mm

Full Set

					
BMT-0010-S With Stainless Steel Housing & Nylon insert	BMT-0020-S With Stainless Steel Housing & Nylon insert	BMT-0030-S With Stainless Steel Housing & Nylon insert	BMT-0040-S With Stainless Steel Housing & Nylon insert	BMT-0050-S With Stainless Steel Housing & Nylon insert	BMT-0060-S With Stainless Steel Housing & Nylon insert

BALL HOUSING MINOR \ BALL NYLON INSERT

		
BHM-5032 Stainless Steel Housing	BHM-4730 Slim Stainless Steel Housing	BIS-0025 Nylon Insert

ANALOG FOR BALL ATTACHMENT


AMB-1201 Analog for Ball attachment

BALL WIDE STRAIGHT







4.5 Platform

					
BWS-4010 H = 1.0 mm Without Cap	BWS-4020 H = 2.0 mm Without Cap	BWS-4030 H = 3.0 mm Without Cap	BWS-4040 H = 4.0 mm Without Cap	BWS-4050 H = 5.0 mm Without Cap	BWS-4060 H = 6.0 mm Without Cap

BALL WIDE SET - NYLON

4.5 Platform

Full Set

					
BWT-0010-S With Stainless Steel Housing & Nylon insert	BWT-0020-S With Stainless Steel Housing & Nylon insert	BWT-0030-S With Stainless Steel Housing & Nylon insert	BWT-0040-S With Stainless Steel Housing & Nylon insert	BWT-0050-S With Stainless Steel Housing & Nylon insert	BWT-0060-S With Stainless Steel Housing & Nylon insert





ANGULAR MULTI SYSTEM CONCEPT 9°, 18°, 30°

Angular Ball Attachment - SET

9°

BALL MINOR ANGULAR Multi System Concept

3.75 Platform

			
BMA-0920-S H02 = 3.2 mm 9° With Stainless Steel Housing & Nylon Insert	BMA-0930-S H03 = 4.2 mm 9° With Stainless Steel Housing & Nylon Insert	BMA-0940-S H04 = 5.2 mm 9° With Stainless Steel Housing & Nylon Insert	BMA-0950-S H05 = 6.2 mm 9° With Stainless Steel Housing & Nylon Insert





ANALOG FOR BALL ATTACHMENT


AMB-1201 Analog for Ball attachment

18°

BALL MINOR ANGULAR Multi System Concept

3.75 Platform

			
BMA-1820-S H02 = 3.5 mm 18° With Stainless Steel Housing & Nylon Insert	BMA-1830-S H03 = 4.5 mm 18° With Stainless Steel Housing & Nylon Insert	BMA-1840-S H04 = 5.5 mm 18° With Stainless Steel Housing & Nylon Insert	BMA-1850-S H05 = 6.5 mm 18° With Stainless Steel Housing & Nylon Insert



BALL HOUSING MINOR \ BALL NYLON INSERT

	
BHM-5032 Stainless Steel Housing	BIS-0025 Nylon Insert
	
BHM-4730 Slim Stainless Steel Housing	

30°

BALL MINOR ANGULAR Multi System Concept

3.75 Platform

			
BMA-3020-S H02 = 4.5 mm 30° With Stainless Steel Housing & Nylon Insert	BMA-3030-S H03 = 5.5 mm 30° With Stainless Steel Housing & Nylon Insert	BMA-3040-S H04 = 6.2 mm 30° With Stainless Steel Housing & Nylon Insert	BMA-3050-S H05 = 7.2 mm 30° With Stainless Steel Housing & Nylon Insert

Designated information for
Angular Multi System Adapters:

Internal Hex

3.75 Platform = EMA-....-M

4.5 Platform = EMA-....-W

TEMPORARY AESTHETIC ABUTMENTS

Includes a Retaining Screw

Abutments for temporary restoration made of biocompatible material called PEEK. Polyether Ether Ether Ketone is a colorless organic ThermoPlastic polymer, proven to provide long term service.

STRAIGHT LEVEL MINOR - PEEK

3.75 Platform

		
SLM-1110-P H = 1.0 mm Shoulder	SLM-1220-P H = 2.0 mm Shoulder	SLM-1330-P H = 3.0 mm Shoulder





STRAIGHT TITANIUM ABUTMENT

3.75 Platform

	
SMT-1210 Engaging	SMT-1210-R Non Engaging

ANGULAR CURVE MINOR - PEEK

3.75 Platform 15°

			
ACM-1115-10-P H = 1.0 mm Aesthetic Shoulder	ACM-1215-20-P H = 2.0 mm Aesthetic Shoulder	ACM-1215-30-P H = 2.5 mm Aesthetic Shoulder	ACM-1315-40-P H = 3.5 mm Aesthetic Shoulder

IMMEDIATE TEMPORARY CONICAL

For immediate loading restoration - For making a temporary bridge

3.75 Platform

						
ITC-0010 H = 1.0 mm Aesthetic Shoulder	ITC-0020 H = 2.0 mm Aesthetic Shoulder	ITC-0030 H = 3.0 mm Aesthetic Shoulder	ITC-0040 H = 4.0 mm Aesthetic Shoulder	ITC-0050 H = 5.0 mm Aesthetic Shoulder	ITC-0060 H = 6.0 mm Aesthetic Shoulder	CCP-00 PEEK cover Non engaging

STRAIGHT AESTHETIC ABUTMENTS

Includes a Retaining Screw





STRAIGHT CURVE MINOR

3.75 Platform

			
SCM-1110 H = 1.0 mm Aesthetic Shoulder	SCM-1220 H = 2.0 mm Aesthetic Shoulder	SCM-1330 H = 3.0 mm Aesthetic Shoulder	SCM-1440 H = 4.0 mm Aesthetic Shoulder

STRAIGHT CURVE CONCAVE MINOR

3.75 Platform

			
SCCM-1110 H = 1.0 mm ConCave Straight Shoulder	SCCM-1220 H = 2.0 mm ConCave Straight Shoulder	SCCM-1330 H = 3.0 mm ConCave Straight Shoulder	SCCM-1440 H = 4.0 mm ConCave Straight Shoulder

STRAIGHT CURVE CONCAVE MINOR **TIN COATING**

3.75 Platform

			
SCCM-1110-T H = 1.0 mm ConCave Straight Shoulder	SCCM-1220-T H = 2.0 mm ConCave Straight Shoulder	SCCM-1330-T H = 3.0 mm ConCave Straight Shoulder	SCCM-1440-T H = 4.0 mm ConCave Straight Shoulder

STRAIGHT LEVEL MINOR WITH HEX

3.75 Platform

		
SLM-1116 H = 1.0 mm Straight Shoulder	SLM-1226 H = 2.0 mm Straight Shoulder	SLM-1336 H = 3.0 mm Straight Shoulder

STRAIGHT LEVEL MINOR

3.75 Platform

		
SLM-1110 H = 1.0 mm Straight Shoulder	SLM-1220 H = 2.0 mm Straight Shoulder	SLM-1330 H = 3.0 mm Straight Shoulder

STRAIGHT LEVEL MINOR WITH GROOVES

3.75 Platform

		
SLM-1110-G H = 1.0 mm Straight Shoulder	SLM-1220-G H = 2.0 mm Straight Shoulder	SLM-1330-G H = 3.0 mm Straight Shoulder




STRAIGHT SLIM MINOR

3.75 Platform

		
SSM-0006 H = 4.0 mm Ultra Slim Abutment	SSM-0008 H = 6.0 mm Ultra Slim Abutment	SSM-0010 H = 8.0 mm Ultra Slim Abutment

STRAIGHT MINOR

3.75 Platform

			
STM-0910 H = 7.0 mm Standard Abutment	STM-1110 H = 9.0 mm Standard Abutment	STM-1310 H = 11.0 mm Standard Abutment	STM-1120 H = 9.0 mm Standard Abutment




STRAIGHT CURVE WIDE

4.5 Platform

			
SCW-1110 h = 1.0 mm Aesthetic Shoulder	SCW-1220 h = 2.0 mm Aesthetic Shoulder	SCW-1330 h = 3.0 mm Aesthetic Shoulder	SCW-1440 h = 4.0 mm Aesthetic Shoulder



STRAIGHT CURVE WIDE

4.5 Platform

			
SCCW-1110 h = 1.0 mm Aesthetic Shoulder	SCCW-1220 h = 2.0 mm Aesthetic Shoulder	SCCW-1330 h = 3.0 mm Aesthetic Shoulder	SCCW-1440 h = 4.0 mm Aesthetic Shoulder

STRAIGHT CURVE WIDE **TIN COATED**

4.5 Platform

			
SCCW-1110-T h = 1.0 mm Aesthetic Shoulder	SCCW-1220-T h = 2.0 mm Aesthetic Shoulder	SCCW-1330-T h = 3.0 mm Aesthetic Shoulder	SCCW-1440-T h = 4.0 mm Aesthetic Shoulder

STRAIGHT LEVEL WIDE

4.5 Platform

		
SLW-1110 h = 1.0 mm Straight Shoulder	SLW-1220 h = 2.0 mm Straight Shoulder	SLW-1330 h = 3.0 mm Straight Shoulder

STRAIGHT WIDE

4.5 Platform


STW-1110 H = 9.0 mm Standard Abutment

H = Abutment height.
h = Shoulder height.

ANGULAR AESTHETIC ABUTMENTS

Includes a Retaining Screw

SPECIFICATIONS

- The angular aesthetic Abutment body functions as impression coping Includes the Prosthetic Screw.
- Installation using large external hex drivers (1.3mm)
- Non-rotational design for standard, internal-hex implants
- Compatible with ZIMMER®, MIS, ALPHA BIO,A.B.DENTAL, Bio Horizon
- Made of Commercial grade 5 Titanium alloy (ASTM-F-136 02A COND A)

In order to preserve the original structure of the titanium alloy, **SURCAM** utilizes advanced CAD/CAM methods in the production process. These methods effectively reduce undesirable deformations of the alloy structure.

ANGULAR CURVE MINOR

3.75 Platform 15°

			
ACM-1115-10 h = 1.0 mm Aesthetic Shoulder	ACM-1215-20 h = 2.0 mm Aesthetic Shoulder	ACM-1215-30 h = 3.0 mm Aesthetic Shoulder	ACM-1315-40 h = 4.0 mm Aesthetic Shoulder




ANGULAR CURVE WIDE

4.50 Platform 15°

	
ACW-1115-10 h = 1.0 mm Aesthetic Shoulder	ACW-1215-20 h = 2.0 mm Aesthetic Shoulder

ANGULAR CURVE CONCAVE MINOR

3.75 Platform 15°

			
ACCM-1115-10 h = 1.0 mm ConCave Shoulder	ACCM-1215-20 h = 2.0 mm ConCave Shoulder	ACCM-1315-30 h = 3.0 mm ConCave Shoulder	ACCM-1415-40 h = 4.0 mm ConCave Shoulder





ANGULAR CURVE CONCAVE TIN COATING

3.75 Platform 15°

			
ACCM-1115-10-T h = 1.0 mm ConCave Shoulder TIN Coating	ACCM-1215-20-T h = 2.0 mm ConCave Shoulder TIN Coating	ACCM-1315-30-T h = 3.0 mm ConCave Shoulder TIN Coating	ACCM-1415-40-T h = 4.0 mm ConCave Shoulder TIN Coating

ANGULAR CURVE MINOR

3.75 Platform 25°

			
ACM-1125-10 h = 1.0 mm Aesthetic Shoulder	ACM-1225-20 h = 2.0 mm Aesthetic Shoulder	ACM-1325-30 h = 3.0 mm Aesthetic Shoulder	ACM-1425-40 h = 4.0 mm Aesthetic Shoulder

H = Abutment height.
h = Shoulder height.

ANGULAR CURVE CONCAVE MINOR

3.75 Platform 25°

			
ACCM-1125-10 h = 1.0 mm ConCave Shoulder	ACCM-1225-20 h = 2.0 mm ConCave Shoulder	ACCM-1325-30 h = 3.0 mm ConCave Shoulder	ACCM-1425-40 h = 4.0 mm ConCave Shoulder

ANGULAR CURVE CONCAVE TIN COATING

3.75 Platform 25°

			
ACCM-1125-10-T h = 1.0 mm ConCave Shoulder TIN Coating	ACCM-1225-20-T h = 2.0 mm ConCave Shoulder TIN Coating	ACCM-1325-30-T h = 3.0 mm ConCave Shoulder TIN Coating	ACCM-1425-40-T h = 4.0 mm ConCave Shoulder TIN Coating

ANGULAR SMOOTH MINOR - LEAF

3.75 Platform 15°

	
ASM-0915-LF H = 7.0 mm Chopped Edge	ASM-1115-LF H = 9.0 mm Chopped Edge

ANGULAR SMOOTH MINOR - LEAF

3.75 Platform 25°

	
ASM-0925-LF H = 7.0 mm Chopped Edge	ASM-1125-LF H = 9.0 mm Chopped Edge

ANGULAR SMOOTH MINOR

3.75 Platform 15°

		
ASM-0915 H = 7.0 mm Smooth Shape	ASM-1115 H = 9.0 mm Smooth Shape	ASM-1315 H = 11.0 mm Smooth Shape

Ultra Thin and LEAF Abutments are for frontal teeth


ANGULAR SMOOTH THIN

3.75 Platform 15° & 25°

	
AST-1115 H = 9.0 mm Ultra Thin 15°	AST-1125 H = 9.0 mm Ultra Thin 25°

ANGULAR SMOOTH MINOR

3.75 Platform 25°


		
ASM-0925 H = 7.0 mm Smooth Shape	ASM-1125 H = 9.0 mm Smooth Shape	ASM-1325 H = 11.0 mm Smooth Shape

ANGULAR SMOOTH WIDE

4.5 Platform 15°


ASW-1115 H = 9.0 mm Smooth Shape

H = Abutment height.
h = Shoulder height.



TRANSFERS, IMPRESSION COPING

Includes a Retaining Screw

TRANSFER MINOR CLOSED
Closed-tray Transfer Abutments

3.75 Platform

				
TMC-0910 Titanium Body H = 1.0 mm	TMC-0920 Titanium Body H = 2.0 mm	TMC-0930 Titanium Body H = 3.0 mm	TMC-0940 Titanium Body H = 4.0 mm	TMC-0950 Titanium Body H = 5.0 mm



TRANSFER MINOR CLOSED
Closed-tray CONCAVE Transfer

3.75 Platform

	
TMCC-0930 Titanium Body CONCAVE H = 3.0 mm	TMCC-0950 Titanium Body CONCAVE H = 5.0 mm


TRANSFER MINOR CLOSED
Closed-tray CONCAVE Transfer

3.75 Platform

	
TMC-100003 Long Transfer L = 6.0 mm Titanium Body	TMC-100004 Long Transfer L = 7.5 mm Titanium Body

TRANSFER MINOR CLOSED
Closed-tray Transfer

3.75 Platform


TMC-0900 Stainless Steel Body

TRANSFER CLOSED-TRAY CLIP

3.75 Platform 4.5 Platform

	
TCMC-1001 Stainless Steel Clip Transfer	TCWC-1001 Stainless Steel Clip Transfer




TRANSFER CLOSED-TRAY CLIP

3.75 Platform 4.5 Platform

	
TCMC-1501 Stainless Steel Clip Transfer	TCWC-1501 Stainless Steel Clip Transfer


TRANSFER MINOR STAINLESS STEEL
Open Tray Transfer

3.75 Platform

		
TMS-1301 Stainless Steel Body Long design	TSS-1001 Stainless Steel Body Slim design	TMS-1001 Stainless Steel Body Short design

TRANSFER WIDE STAINLESS STEEL
Open Tray Transfer

4.5 Platform


TWS-1301 Stainless Steel Body Long design

TRANSFER WIDE CLOSED
Closed-tray Transfer Abutments

4.5 Platform

				
TWC-0910 Titanium Body H = 1.0 mm	TWC-0920 Titanium Body H = 2.0 mm	TWC-0930 Titanium Body H = 3.0 mm	TWC-0940 Titanium Body H = 4.0 mm	TWC-0950 Titanium Body H = 5.0 mm

TRANSFER MATCHED PLASTIC
Closed-tray Transfer Plastic Cap

		
TMP-0008 Low Level Plastic Cap	TMP-0010 High Level Plastic Cap	CPC-0008 Burn Out Plastic Cap

PLASTIC ABUTMENTS FOR CASTING

STRAIGHT PLASTIC FOR CASTING

3.75 Platform \ 4.5 Platform

			
PSM-0013 Round Base	PSW-0013 Round Base Wide Platform	PSM-6013 Hex Base	PSW-6013 Hex Base Wide Platform

TITANIUM PLASTIC ABUTMENT

3.75 Platform

	
TPA-1420 Hex Base	TPA-1420-R Round Base

TITANIUM PLASTIC ABUTMENT

4.5 Platform

	
TPA-1420-W Hex Base Wide Platform	TPA-1420-WR Round Base Wide Platform

CrCo PLASTIC ABUTMENT

3.75 Platform

	
CPA-1405 Hex Base H= 0.5 mm	CPA-1405-R Round Base H= 0.5 mm

CrCo PLASTIC ABUTMENT

3.75 Platform

	
CPA-1420 Hex Base H= 2.0 mm	CPA-1420-R Round Base H= 2.0 mm

CrCo PLASTIC ABUTMENT

3.75 Platform

	
CPA-1430 Hex Base H= 3.0 mm	CPA-1430-R Round Base H= 3.0 mm

CrCo PLASTIC ABUTMENT


4.5 Platform

	
CPA-1420-W Hex Base Wide Platform	CPA-1420-WR Round Base Wide Platform

The CrCo / Titanium / Plastic Abutments are for investment and casting following conventional crown and bridge techniques. The CrCo / Titanium / Plastic Abutments combine the precision of a machined interface between the implant and the Abutment with the convenience of a Plastic sleeve for casting.

Analogs are used for the fabrication of a master model. They are replicas of the implant. All analogs are made of high grade stainless steel.

ANALOG STAINLESS STEEL



ANALOG FOR MULTI-UNIT



ANALOG FOR WISE CLICK SYSTEM



ANALOG FOR BALL ATTACHMENT



3D SCANNING TiBase

3.75 Platform



SMS-0040
CadCam TiBase
Rotational



SMS-6040
CadCam TiBase
Non-rotational

4.5 Platform



SWS-0040
CadCam TiBase
wide
Rotational



SWS-6040
CadCam TiBase
wide
Non-rotational



SBT-100016
3D Scanbody
for TiBase / PEEK
Non-rotational



SCSC-1720
Long Screw for
3D Scanning

AMANN GIRRBACH®



BMZ-6001
CadCam TiBase
Non-rotational



BMZ-0001
CadCam TiBase
Rotational

SIRONA®



SMS-6050
CadCam TiBase
sro. Non-rotational



SMS-6051
CadCam TiBase
sro. Non-rotational



SMS-6052
CadCam TiBase
sro. Non-rotational



SMS-6053
CadCam TiBase
sro. Non-rotational



SMS-6054
CadCam TiBase
sro. Non-rotational

SIRONA®



SMS-6050-R
CadCam TiBase
sro. Rotational

ZIRKONZAHN®



SMS-0030
CadCam TiBase zrż
Rotational



SMS-6030
CadCam TiBase zrż
Non-rotational

ZIRKONZAHN®



SWS-0030
CadCam TiBase
wide
Rotational



SWS-6030
CadCam TiBase
wide
Non-rotational

MIS® on ZIRKONZAHN®



SNS-400002
CadCam TiBase NP
Rotational



SNS-400004
CadCam TiBase NP
Non-rotational

MIS® on SIRONA®



SNS-400007
CadCam TiBase NP
Non-rotational



SNS-400008
CadCam TiBase NP
Rotational

MIS®




SMS-100001
CadCam TiBase
Non-rotational



SMS-100002
CadCam TiBase
Rotational

INTERNAL EXTERNAL ADAPTER

External Hex (Top): Regular 4.1
Internal Hex (Bottom): Regular 3.75





	
INEX-2427-R Round Base 3.75 Platform	INEX-2427 Hex Base 3.75 Platform

Related Screws

SCIE-2510	SCIE-2510
-----------	-----------

INTERNAL EXTERNAL ADAPTER

External Hex (Top): Wide 5.0 Internal Hex (Bottom): Wide 4.5	External Hex (Top): Regular 4.1 Internal Hex (Bottom): Wide 4.5
---	--

			
INEX-2433-WR Round Base 4.5 Platform	INEX-2433-W Hex Base 4.5 Platform	INEX-2427-WR Round Base 4.5 Platform	INEX-2427-W Hex Base 4.5 Platform







Related Screws

SCIE-3009	SCIE-3009	SCIE-2510	SCIE-2510
-----------	-----------	-----------	-----------

Healing caps are used to maintain access through tissue during the healing phase. They are available in a variety of heights and diameters to accommodate specific tissue depth.







HEALING CAP MINOR ULTRA SLIM

3.75 Platform

					
HCM-3720 Dia. 3.7 mm H = 2.0 mm	HCM-3730 Dia. 3.7 mm H = 3.0 mm	HCM-3740 Dia. 3.7 mm H = 4.0 mm	HCM-3750 Dia. 3.7 mm H = 5.0 mm	HCM-3760 Dia. 3.7 mm H = 6.0 mm	HCM-3770 Dia. 3.7 mm H = 7.0 mm







HEALING CAP MINOR SLIM

3.75 Platform

					
HCM-4020 Dia. 4.6 mm H = 2.0 mm	HCM-4030 Dia. 4.6 mm H = 3.0 mm	HCM-4040 Dia. 4.6 mm H = 4.0 mm	HCM-4050 Dia. 4.6 mm H = 5.0 mm	HCM-4060 Dia. 4.6 mm H = 6.0 mm	HCM-4070 Dia. 4.6 mm H = 7.0 mm

HEALING CAP MINOR NORMAL

3.75 Platform

					
HCM-5020 Dia. 5.0 mm H = 2.0 mm	HCM-5030 Dia. 5.0 mm H = 3.0 mm	HCM-5040 Dia. 5.0 mm H = 4.0 mm	HCM-5050 Dia. 5.0 mm H = 5.0 mm	HCM-5060 Dia. 5.0 mm H = 6.0 mm	HCM-5070 Dia. 5.0 mm H = 7.0 mm

CONCAVE HEALING CAP

3.75 Platform

			
CHC-0155 Dia. 5.5 mm H = 1.0 mm	CHC-0255 Dia. 5.5 mm H = 2.0 mm	CHC-0355 Dia. 5.5 mm H = 3.0 mm	CHC-0455 Dia. 5.5 mm H = 4.0 mm

HEALING CAP WIDE

4.5 Wide Platform

					
HCW-6030 Dia. 6.0 mm H = 3.0 mm	HCW-6040 Dia. 6.0 mm H = 4.0 mm	HCW-6050 Dia. 6.0 mm H = 5.0 mm	HCW-6060 Dia. 6.0 mm H = 6.0 mm	HCW-6070 Dia. 6.0 mm H = 7.0 mm	HCW-6080 Dia. 6.0 mm H = 8.0 mm

SUPREME HEALING CAP

3.75 Platform

		
HCS-37-- Dia. 3.7 mm H = 2.0-7.0 mm	HCS-40-- Dia. 4.6 mm H = 2.0-7.0 mm	HCS-50-- Dia. 5.0 mm H = 2.0-7.0 mm

Supreme Healing Cap With TiN Coating




The TiN coatings modifies only surface properties and has no effect on the substrate properties or biomechanical functionality.

Biocompatibility:

Biocompatibility of TiN has been demonstrated in numerous published studies. It is known to Improve friction and has a segnificantly higher degree of surface scratch resistance compared to other materials.

PROSTHETIC SCREWS


SCREWS FOR ABUTMENTS

		
SCA-2385 For Angular Abutment	SCS-3095 For Straight Abutment	SMU-4023 For Multi Unit Sleeve

SCREWS FOR OVERDENTURE ABUTMENTS

		
SCS-3011 For Overdenture Abutment ECA-6005	SCS-3012 For Overdenture Abutment ECA-6015	SCS-3013 For Overdenture Abutment ECA-6025

SCREWS FOR OPEN TRAY TRANSFER

		
SCT-3018 Screw for Short Open tray transfer	SCT-3023 Screw for Long open tray transfer	SCSC-1720 Long Screw for 3D scanning

SCREW FOR Plastic ABUTMENTS

	
SCP-2580 For Straight Plastic Abutment	SCA-2565 For Angular Plastic Abutment

WRENCHES

&

COMPONENTS





Manual Tools



Machine Tools



Ratchet Hex
6.35 Wrench



Corckskey



Adapter

96 ◀ Manual Tools

Machine Tools ▶ 97

98 ◀ Ratchet Hex 6.35 Wrench

99 ◀ Corckskey

Adapter ▶ 100

MANUAL TOOLS


PROSTHETIC SCREW DRIVES 1.27

Manual

			
KSH-1013 10.0 mm Projection	KSH-1513 15.0 mm Projection	KSH-2013 20.0 mm Projection	KSH-3013 30.0 mm Projection

INTERNAL HEX IMPLANT KEY WITH HOLDING CLIP

Manual

	
KSH-1224 12.0 mm Projection	KSH-1824 18.0 mm Projection

PROSTHETIC SCREW DRIVES 1.27

Machine

		
MKS-0813 8.0 mm Projection	MKM-1513 15.0 mm Projection	MKL-2513 25.0 mm Projection

INTERNAL HEX IMPLANT KEY WITH HOLDING CLIP


Machine

	
MIC-1224 12.0 mm Projection	MIC-1824 18.0 mm Projection


RATCHET WRENCH

REGULAR RATCHET WRENCH HEX DRIVE 6.35 mm

It is recommended to use over a Torque Wrench



Tightening with the correct torque minimizes the risk of the screw loosening and improves the precision during bridge fabrication. Torque control can also reduce the risk of stress being built into the restoration and prevents component fractures.



CORCKSKEY

CORCKSKEY


Universal Screw and Implant Extractor



Can be used to remove:

- Implants with external hex connection , Internal three leg fit ("replace select") connection and Internal hex connection.
- Prosthetic screw with damage driver.

SCREW RELEASE KEY



THE CORCKSKEY HASE



ADAPTER

SUPREME ADAPTER FOR HANDLE DRIVERS



AHD-10
Standard Tightening
Adapter. 10.0 mm

ADAPTER FOR HANDLE DRIVERS
Fastening strength increases



AHD-12
Extra increased tightening
Adapter. 12.0 mm

PRO
POSITIONS



The U-Type


Classic Thread Conical shape.

Internal Hex 11° Conical Fit Platform.

The U-Type Implant is suitable for all clinical procedures.

With the U-Type Implant you achieve high primary stability.

Moreover, in cases of weak bone structure the U-Type Implant condenses the peri-implant bone. The connection platform has high mechanical strength and a double hex that offers 12 prosthetic positions with a 30° deviation. The U-Type Implant is very stable during placement with the option of immediate restoration and loading.



	A	B
U35008CF	3.50	2.20
U350100CF	3.50	2.20
U350115CF	3.50	2.20
U350013CF	3.50	2.20


	A	B
U40008CF	4.00	2.40
U400100CF	4.00	2.40
U400115CF	4.00	2.40
U40013CF	4.00	2.40

The C-Type

Classic Thread Parallel shape.

Internal Hex 11° Conical Fit Platform.

The C-Type Implant's cylindrical body is designed to maximize initial stability in soft bone. The shallow cutting edges on its threads are specifically designed for soft bone indications and minimize pressure on hard bone. The connection platform has high mechanical strength and a double hex that offers 12 prosthetic positions with a 30° deviation. The C-Type Implant is very stable during placement with the option of immediate restoration and loading.




	A	B
C35008CF	3.50	2.20
C350100CF	3.50	2.20
C350115CF	3.50	2.20
C35013CF	3.50	2.20

	A	B
C40008CF	4.00	2.40
C400100CF	4.00	2.40
C400115CF	4.00	2.40
C40013CF	4.00	2.40

	A	B
C45008CF	4.50	2.91
C450100CF	4.50	2.91
C450115CF	4.50	2.91
C45013CF	4.50	2.91

	A	B
C50008CF	5.00	3.41
C500100CF	5.00	3.41
C500115CF	5.00	3.41
C50013CF	5.00	3.41

Surcam's Implants are packed in sterile tubes where they remain until required. The implants are protected by a specially designed stainless steel "Implant carrier" which both protects them physically, and is also a tool that guides the dental surgeon to insert them in the patient's mouth at the desired orientation. The implant carrier support max. 40Ncm torque and should be replaced with a suitable tool for higher torques if required. This is done without any direct contact with the implant surface.



U-TYPE CONICAL FIT

CONICAL FIT PLATFORM 3.5 / 4.0

Dia. 3.5 mm

			
U35008CF L = 8 mm	U350100CF L = 10 mm	U350115CF L = 11.5 mm	U350013CF L = 13 mm

Dia. 4.0 mm

			
U40008CF L = 8 mm	U400100CF L = 10 mm	U400115CF L = 11.5 mm	U40013CF L = 13.0 mm


C-TYPE CONICAL FIT

CONICAL FIT PLATFORM 3.5 / 4.0

Dia. 3.5 mm





			
C35008CF L = 8 mm	C350100CF L = 10 mm	C350115CF L = 11.5 mm	C350013CF L = 13 mm

COVER SCREW


SCC-3540 Cover Screw 3.5 / 4.0 Platform

CONICAL FIT PLATFORM 3.5 / 4.0

Dia. 4.0 mm

			
C40008CF L = 8 mm	C400100CF L = 10 mm	C400115CF L = 11.5 mm	C40013CF L = 13 mm

CONICAL FIT PLATFORM 3.5 / 4.0

Dia. 4.5 mm


			
C45008CF L = 8 mm	C450100CF L = 10 mm	C450115CF L = 11.5 mm	C45013CF L = 13 mm

CONICAL FIT PLATFORM 3.5 / 4.0

Dia. 5.0 mm

			
C50008CF L = 8 mm	C500100CF L = 10 mm	C500115CF L = 11.5 mm	C50013CF L = 13 mm

Drilling Sequence - U-Type Conical fit




● Drill to part depth of the Implant (*See illustrations 1 & 2)

○ Drill to full depth of the Implant (*See illustrations 1 & 2)

▶ Drill only through the cortical bone, the drill should not be used to full depth.


The chart above illustrates the recommended drilling sequence for SURCAM's implants. Clinical judgement with respect to each individual patient must supersede any recommendation made.

illustrations 1




For bone type I II

illustrations 2




For bone type III IV

Drilling Sequence - C-Type Conical fit




- Drill to part depth of the Implant (*See illustrations 1 & 2)
- Drill to full depth of the Implant (*See illustrations 1 & 2)

 Drill only through the cortical bone, the drill should not be used to full depth.
Keep 2.5mm of the bottom hole with the 2.0mm Dia. drill.


The chart above illustrates the recommended drilling sequence for SURCAM's implants. Clinical judgement with respect to each individual patient must supersede any recommendation made.

illustrations 1



For **bone type I II**

illustrations 2



For **bone type III IV**

2.50 mm
Small Dia. Drill

CONICAL FIT HEALING CAP

Platform 3.5 / 4.0 mm

					
HCC-4502 Dia. 4.5 mm H = 2.0 mm	HCC-4504 Dia. 4.5 mm H = 4.0 mm	HCC-4506 Dia. 4.5 mm H = 6.0 mm	HCC-5502 Dia. 5.5 mm H = 2.0 mm	HCC-5504 Dia. 5.5 mm H = 4.0 mm	HCC-5506 Dia. 5.5 mm H = 6.0 mm

CONICAL FIT TRANSFER

Closed tray



	
TCF-0930 Titanium Abutment L = 3.0 mm	TCF-0950 Titanium Abutment L = 5.0 mm

TRANSFER MATCHED Plastic
Closed-tray Transfer Plastic Cap

		
TMP-0008 Low Level Plastic Cap	TMP-0010 High Level Plastic Cap	CPC-0008 Burn Out Plastic Cap



Open tray

	
TSF-1405 Stainless-steel Abutment L = 14.0 mm	TLF-1905 Stainless-steel Abutment L = 19.0 mm

ANALOG




ACS-3540 3.5 / 4.0 Platform

CONICAL FIT PROSTHETIC ABUTMENT

STRAIGHT LEVEL

	
SLS-1115 h = 1.5 mm	SLS-1330 h = 3.0 mm
REF ASTRA® 24285	REF ASTRA® 24286

15° ANGULAR LEVEL

	
ALS-4515-35 h = 3.5 mm	ALS-5515-35 h = 3.5 mm
REF ASTRA® 24291	REF ASTRA® 24292

* The original ASTRA abutment is with 20° angulation.




CHROME COBALTE ABUTMENT FOR CASTING

CrCo PLASTIC ABUTMENT




3.5 / 4.0 Platform

		
ACC-200001 CrCo base RP hex h = 1 mm	ACC-200002 CrCo base RP hex h = 2 mm	ACC-200003 CrCo base RP hex h = 3 mm

CONICAL FIT BALL ATTACHMENT

		
BCF-0010 L=1.0 mm	BCF-0020 L=2.0 mm	BCF-0030 L=3.0 mm

SET

		
BCF-0010-S L=1.0 mm	BCF-0020-S L=2.0 mm	BCF-0030-S L=3.0 mm

CONICAL FIT IMPLANT KEYS

HAND TOOL


		
HIS-1214 Hand CF implant key - SHORT L=14 mm	HIM-1218 Hand CF implant key MEDIUM L=18 mm	HIL-1222 Hand CF implant key LONG L=22 mm

MACHINE TOOL

	
MIS-1209 Machine CF implant key - SHORT L=9 mm	MIM-1214 Machine CF implant key - MEDIUM L=14 mm

INTERNAL

THREE LEG FIT





Compatible - Internal Three Leg Fit platform: 3.5, 4.3 & 5.0 mm Implants.
[Noble Biocare®, Replace Select®, Tri-lob®].


INVERSE ANGULAR ABUTMENTS

Introduction & Advantages


The **Inverse** Abutments allow 6 different positioning options for better flexibility of implant positioning.



120° deviation
allows only
3 different positions



By using 2 Abutments (Regular and Inverse) we achieve a 60° deviation which allows 6 different positioning options for better flexibility of implant positioning.



Components
for
OverDenture
Restoration



Aesthetic
Abutments

Components
for
Cement
Restoration



Casting
Abutments



Analog




3D Scanning
Interface

Components
for
Lab
Use



Healing
Cap

Complementary
Items



STRAIGHT WISE CLICK

Overdenture Attachment connected to Implants

WISE CLICK REGULAR - RETAINED OVER DENTURE ABUTMENTS
Abutment bases connected to implants

4.3 Regular Platform

For Full Set add SPE to catalog number

					
WCR-0000 H = 0.3 mm	WCR-0010 H = 1.0 mm	WCR-0020 H = 2.0 mm	WCR-0030 H = 3.0 mm	WCR-0040 H = 4.0 mm	WCR-0050 H = 5.0 mm

WISE CLICK REGULAR - SET PEEK

Abutment bases with peek attachment and housing covers

4.3 Platform

FULL SET

					
WCR-0000-SPE H = 0.3 mm	WCR-0010-SPE H = 1.0 mm	WCR-0020-SPE H = 2.0 mm	WCR-0030-SPE H = 3.0 mm	WCR-0040-SPE H = 4.0 mm	WCR-0050-SPE H = 5.0 mm

WISE CLICK - ACCESSORIES

Included in SET

						
WCH-5628 Stainless Steel Housing Cover	WCH-4823 Stainless Steel Housing	WCP-2348-H Hard Housing PEEK Attachment	WCPF-2348-H Positioning PEEK Attachment	WCPS-2348-H Soft Housing PEEK Attachment	WCP-2348 Direct PEEK Attachment	SCH-8540 Nylon Sealed Ring

Not Included in SET

		
TWC-1105 Transfer For Wise Click	WCT-900 Wise Click Tool	AWC-1301 Analog For Wise Click

ANGULAR WISE CLICK CONCEPT 9°, 18°, 30°

Wise Click Overdenture Attachment - SETS

9°

WISE CLICK ANGULAR - MULTI SYSTEM CONCEPT

Wise Click angular over denture attachment - Set

			
WCA-0920-R H02 = 3.2 mm 9° Wise Click Angular Multi System	WCA-0930-R H03 = 4.2 mm 9° Wise Click Angular Multi System	WCA-0940-R H04 = 5.2 mm 9° Wise Click Angular Multi System	WCA-0950-R H05 = 6.2 mm 9° Wise Click Angular Multi System

18°

WISE CLICK ANGULAR - MULTI SYSTEM CONCEPT





Wise Click angular overdenture Attachment - Set

			
WCA-1820-R H03 = 3.5 mm 18° Wise Click Angular Multi System	WCA-1830-R H04 = 4.5 mm 18° Wise Click Angular Multi System	WCA-1840-R H05 = 5.5 mm 18° Wise Click Angular Multi System	WCA-1850-R H06 = 6.5 mm 18° Wise Click Angular Multi System

30°


WISE CLICK ANGULAR - MULTI SYSTEM CONCEPT

Wise Click angular overdenture Attachment - Set

			
WCA-3020-R H03 = 4.5 mm 30° Wise Click Angular Multi System	WCA-3030-R H04 = 5.5 mm 30° Wise Click Angular Multi System	WCA-3040-R H05 = 6.2 mm 30° Wise Click Angular Multi System	WCA-3050-R H06 = 7.2 mm 30° Wise Click Angular Multi System

MULTI SYSTEM CONCEPT


Angular Correction System



↓
Multi-Unit




↓
Ball Attachment



Note

Full description and information about the accessories find in pages # 72+77.


↓
Wise Click



The **MULTI SYSTEM** is designed to fit all angular over denture applications. These unique Multi System Adapters extend the application range of SURCAM's overdenture lines of Angular Ball Attachment; Angular Multi-Unit and Angular Wise Click System.



Dentures with
ANGULAR MULTI-UNIT



Removable Dentures with
ANGULAR BALL ATTACHMENT



Removable Dentures with
AN ANGULAR WISE CLICK SYSTEM







COMPATIBLE AESTHETIC ABUTMENTS

Internal Three Leg Fit

AESTHETIC REGULAR THREE UNIT - SYSTEM ADAPTER

Angular Multi System Adapter 9°, 18°, 30°

4.3 Regular Platform

INVERSE					
					
ERTU-0943 9° Angle Multi System Adapter	ERTU-1843 18° Angle Multi System Adapter	ERTU-3043 30° Angle Multi System Adapter	ERTU-0943-I 9° Angle Multi System Adapter	ERTU-1843-I 18° Angle Multi System Adapter	ERTU-3043-I 30° Angle Multi System Adapter

WISE CLICK CONNECTION

Wise Click connection for multi system concept

			
WCC-0020 2 mm Connection Multi System Wise Click	WCC-0030 3 mm Connection Multi System Wise Click	WCC-0040 4 mm Connection Multi System Wise Click	WCC-0050 5 mm Connection Multi System Wise Click

WISE CLICK - ACCESSORIES

Included in SET





						
WCH-5628 Stainless Steel Housing Cover	WCH-4823 Stainless Steel Housing	WCP-2348-H Hard Housing PEEK Attachment	WCPF-2348-H Positioning PEEK Attachment	WCPS-2348-H Soft Housing PEEK Attachment	WCP-2348 Direct PEEK Attachment	SCH-8540 Nylon Sealed Ring

Not Included in SET

		
TWC-1105 Transfer For Wise Click	WCT-900 Wise Click Tool	AWC-1301 Analog For Wise Click

BALL MULTI ANGULAR

Ball Attachment connection for multi system concept

			
BMA-0020 2 mm Connection Multi System Ball	BMA-0030 3 mm Connection Multi System Ball	BMA-0040 4 mm Connection Multi System Ball	BMA-0050 5 mm Connection Multi System Ball

BALL HOUSING MINOR \ BALL NYLON INSERT



		
BHM-5032 Stainless Steel Housing	BHM-4730 Slim Stainless Steel Housing	BIS-0025 Nylon Insert

ANALOG FOR BALL ATTACHMENT


AMB-1201 Analog for Ball attachment

MULTI-UNIT CONNECTION

Screw-in connection for multi system concept

			
MUC-0021 2 mm Connection Multi System	MUC-0031 3 mm Connection Multi System	MUC-0041 4 mm Connection Multi System	MUC-0051 5 mm Connection Multi System




MULTI-UNIT STRAIGHT AESTHETIC ABUTMENTS



Abutments with smooth, non-locking surface for Multi-Unit Restoration

AESTHETIC REGULAR UNIT

Straight Multi-Unit

4.3 Regular Platform

		
SRTU-4310 H = 1.0 mm	SRTU-4320 H = 2.0 mm	SRTU-4330 H = 3.0 mm



MULTI SYSTEM



See the instructions Manual

NOTE: The screw geometry was carefully designed and optimized to allow safe tightening and long service. Special attention must be paid while tightening the screw. Make sure to use an original SURCAM hex key with edges that do not show wear. It is recommended to purchase additional screws

MULTI-UNIT - ACCESSORIES

							
ISU-1201 H01	ISU-1202 H02	ISU-1203 H03	ISU-1801 H01	ISU-1802 H02	ISU-1803 H03	IBB-2530 Ball attachment	IBM-3230 Connection bar to be modified to the space between the sleeves
Ti Sleeve with 3 positioning options			Ti Sleeve with 2 positioning options				

							
PSU-0032 Plastic Sleeve screw included	TSU-1250 Titanium Sleeve screw included L = 12.0 mm	CMU-100006 Short slim titanium sleeve	CMU-100007 Long slim titanium sleeve L = 8.0 mm	PSU-1050 3D Scanning Lab device	PSU-100015 3D Scanning Intraoral device	PSU-0050 Parallel connection	TSM-100013 Titanium sleeve for welding

						
TMU-0532 Titanium Short Sleeve screw included	TMU-1000 Open Tray Short Transfer	TMU-1200 Open Tray Transfer screw included	TMU-8005 Closed-Tray Transfer	WQP-011 Plastic for Closed-Tray Transfer	HCC-4550 Healing Cap for Multi-Unit	AMU-1201 Analog for Multi-Unit

ANGULAR MULTI SYSTEM CONCEPT 9°, 18°, 30°





Multi-Unit Overdenture Attachment

9°

ANGULAR MULTI-UNIT SYSTEM CONCEPT

Aesthetic Multi-Unit Angular Overdenture Attachment

4.3 Regular Platform

			
EMA-0921-R H02 = 3.2 mm 9° Angular Multi System	EMA-0931-R H03 = 4.2 mm 9° Angular Multi System	EMA-0941-R H04 = 5.2 mm 9° Angular Multi System	EMA-0951-R H05 = 6.2 mm 9° Angular Multi System

18°

ANGULAR MULTI-UNIT SYSTEM CONCEPT

Aesthetic Multi-Unit Angular Overdenture Attachment

4.3 Regular Platform

			
EMA-1831-R H03 = 3.5 mm 18° Angular Multi System	EMA-1841-R H04 = 4.5 mm 18° Angular Multi System	EMA-1851-R H05 = 5.5 mm 18° Angular Multi System	EMA-1861-R H06 = 6.5 mm 18° Angular Multi System

30°

ANGULAR MULTI-UNIT SYSTEM CONCEPT

Aesthetic Multi-Unit Angular Overdenture Attachment

4.3 Regular Platform

			
EMA-3031-R H03 = 4.5 mm 30° Angular Multi System	EMA-3041-R H04 = 5.5 mm 30° Angular Multi System	EMA-3051-R H05 = 6.2 mm 30° Angular Multi System	EMA-3061-R H06 = 7.2 mm 30° Angular Multi System

ANGULAR BALL ATTACHMENT - SET

Overdenture Ball Attachment connection to Implants

9°

BALL REGULAR ANGULAR MULTI SYSTEM CONCEPT





4.3 Regular Platform

			
BRA-0920-S H02 = 3.2 mm With Stainless Steel Housing & Silicon Insert	BRA-0930-S H03 = 4.2 mm With Stainless Steel Housing & Silicon Insert	BRA-0940-S H04 = 5.2 mm With Stainless Steel Housing & Silicon Insert	BRA-0950-S H05 = 6.2 mm With Stainless Steel Housing & Silicon Insert

18°

BALL REGULAR ANGULAR MULTI SYSTEM CONCEPT





4.3 Regular Platform

			
BRA-1820-S H03 = 3.5 mm With Stainless Steel Housing & Silicon Insert	BRA-1830-S H04 = 4.5 mm With Stainless Steel Housing & Silicon Insert	BRA-1840-S H05 = 5.5 mm With Stainless Steel Housing & Silicon Insert	BRA-1850-S H06 = 6.5 mm With Stainless Steel Housing & Silicon Insert

30°

BALL REGULAR ANGULAR MULTI SYSTEM CONCEPT

4.3 Regular Platform

			
BRA-3020-S H03 = 4.5 mm With Stainless Steel Housing & Silicon Insert	BRA-3030-S H04 = 5.5 mm With Stainless Steel Housing & Silicon Insert	BRA-3040-S H05 = 6.2 mm With Stainless Steel Housing & Silicon Insert	BRA-3050-S H06 = 7.2 mm With Stainless Steel Housing & Silicon Insert

BALL REGULAR STRAIGHT

4.3 Regular Platform

		
BRST-4010 H = 1.0 mm Without Cap	BRST-4020 H = 2.0 mm Without Cap	BRST-4030 H = 3.0 mm Without Cap

BALL REGULAR SET - SILICON

4.3 Regular Platform

		
BRST-0010-S Without Stainless Steel Housing & Silicon Insert	BRST-0020-S Without Stainless Steel Housing & Silicon Insert	BRST-0030-S Without Stainless Steel Housing & Silicon Insert

BALL HOUSING MINOR \
BALL NYLON INSERT

		
BHM-5032 Stainless Steel Housing	BHM-4730 Slim Stainless Steel Housing	BIS-0025 Nylon Insert

ANALOG FOR BALL ATTACHMENT


AMB-1201 Analog for Ball attachment

AESTHETIC ABUTMENTS

STRAIGHT CURVE NARROW THREE

3.5 Narrow Platform

	
SCNT-3500-10 1.0 mm Aesthetic Shoulder	SCNT-3500-30 3.0 mm Aesthetic Shoulder

STRAIGHT CURVE REGULAR THREE

4.3 Regular Platform

	
SCRT-4300-10 1.0 mm Aesthetic Shoulder	SCRT-4300-30 3.0 mm Aesthetic Shoulder

STRAIGHT CURVE WIDE THREE



5.0 Wide Platform

	
SCWT-5000-10 1.0 mm Aesthetic Shoulder	SCWT-5000-25 2.5 mm Aesthetic Shoulder

AESTHETIC ABUTMENTS

ANGULAR CURVE NARROW THREE

3.5 Narrow Platform

	
ACNT-3515-10 1.0 mm Aesthetic Shoulder	ACNT-3515-30 3.0 mm Aesthetic Shoulder

INVERSE	
	
ACNT-3515-10-I 1.0 mm Aesthetic Shoulder	ACNT-3515-30-I 3.0 mm Aesthetic Shoulder

ANGULAR CURVE REGULAR THREE

4.3 Regular Platform

	
ACRT-4315-10 1.0 mm Aesthetic Shoulder	ACRT-4315-30 3.0 mm Aesthetic Shoulder

INVERSE	
	
ACRT-4315-10-I 1.0 mm Aesthetic Shoulder	ACRT-4315-30-I 3.0 mm Aesthetic Shoulder

ANGULAR CURVE WIDE THREE

5.0 Wide Platform

	
ACWT-5015-05 1.0 mm Aesthetic Shoulder	ACWT-5015-30 3.0 mm Aesthetic Shoulder

INVERSE	
	
ACWT-5015-05-I 1.0 mm Aesthetic Shoulder	ACWT-5015-30-I 3.0 mm Aesthetic Shoulder

The Titanium / Plastic Abutment are for Investment and casting following conventional crown and bridge techniques. The Titanium / Plastic Abutments combine the precision of a machined interface between the implant and the abutment, with the convenience of a plastic sleeve for casting.

TITANIUM PLASTIC ABUTMENT

4.3 Regular Platform




	
TPAR-1420 Hex Base Regular Platform	TPAR-1420-R Round Base Regular Platform

5.0 Wide Platform

	
TPAR-1420-W Hex Base Wide Platform	TPAR-1420-WR Round Base Wide Platform

Analogs are used for the fabrication of a master model.
They are replicas of the implant. All analogs are made of high grade stainless steel.

ANALOG COMPATIBLE INTERNAL THREE Leg Fit

		
ANS-1135 3,5 mm Narrow Platform	ARS-1143 4,3 mm Regular Platform	AWS-1150 5,0 mm Wide Platform

ANALOG FOR MULTI-UNIT


AMU-1201 Analog for Multi-Unit

ANALOG FOR WISE CLICK SYSTEM


AWC-1301 Analog for Wise Click

ANALOG FOR BALL ATTACHMENT


AMB-1201 Analog for Ball attachment

3D SCANNING INTERFACE

STRAIGHT SCANNING ABUTMENTS

Scanning Titanium Abutments for 3D
Suitable as a Zirconium Interface with the Implant

		
SNS-3530 3.5 mm Narrow Platform	SRS-4330 4.3 mm Regular Platform	SWS-5030 5.0 mm Wide Platform




3.75 Platform

			
SMS-6050 Hex Base SIRONA®	SMS-0050 Round Base SIRONA®	BMZ-6001 Hex Base For Zirconium Abutment	BMZ-0001 Round Base For Zirconium Abutment

INTERNAL EXTERNAL ADAPTER

4.3 Regular Platform

	
INEX-2743-R Round Base	INEX-2743 Internal Three Leg Fit

INTERNAL EXTERNAL ADAPTER

5.0 Wide Platform

	
INEX-2750-R Round Base	INEX-2750 Internal Three Leg Fit

SCREW FOR INTERNAL / EXTERNAL

4.3 Regular Platform



5.0 Wide Platform



Healing caps are used to maintain access through tissue during the healing phase. They are available in a variety of heights and diameters to accommodate specific tissue depth.

WIDE HEALING CAP REGULAR

4.3 Wide Platform / 5.3 mm Dia.


	
HCRT-5330 Dia. = 5.3 mm H = 3.0 mm	HCRT-5350 Dia. = 5.3 mm H = 5.0 mm
Recommended for Angular Restoration	

SLIM HEALING CAP REGULAR


4.3 Regular Platform / 4.3 mm Dia.

	
HCRT-4330 Dia. = 4.3 mm H = 3.0 mm	HCRT-4350 Dia. = 4.3 mm H = 5.0 mm
Recommended for Narrow Spaces	

EXTERNAL HEX





Compatible with external hex platform: 4.1 & 5.0 mm Implants.
[Brenemark®].




Components
for
Overdenture
Restoration



Components
for
Cement
Restoration



Components
for
Lab
Use



Complementary
Items

140 ◀ Wise Click

Multi-Unit
Restoration ▶ 141

142 ◀ Ball
Attachment

143 ◀ Abutments

Casting
Abutments ▶ 144

145 ◀ Analog

Healing
Cap ▶ 146

147 ◀ Transfers

STRAIGHT WISE CLICK

Overdenture Attachment connection to Implants

WISE CLICK FOR EXTERNAL HEX IMPLANTS

Abutment bases connection to implants

4.1 Regular Platform

For Full Set add SPE to catalog number

					
WCRE-4110 H = 1.5 mm	WCRE-4120 H = 2.5 mm	WCRE-4130 H = 3.5 mm	WCRE-4140 H = 4.5 mm	WCRE-4150 H = 5.5 mm	WCRE-4160 H = 6.5 mm

Recommended Tightening Torque 20-25 NCM

WISE CLICK REGULAR - SET PEEK

Abutment bases with PEEK attachment and housing covers

4.1 Platform

FULL SET

					
WCRE-4110-SPE H = 1.0 mm	WCRE-4120-SPE H = 2.0 mm	WCRE-4130-SPE H = 3.0 mm	WCRE-4140-SPE H = 4.0 mm	WCRE-4150-SPE H = 5.0 mm	WCRE-4160-SPE H = 6.0 mm

WISE CLICK - ACCESSORIES

Included in SET

						
WCH-5628 Stainless Steel Housing Cover	WCH-4823 Stainless Steel Housing	WCP-2348-H Hard Housing PEEK Attachment	WCPF-2348-H Positioning PEEK Attachment	WCPS-2348-H Soft Housing PEEK Attachment	WCP-2348 Direct PEEK Attachment	SCH-8540 Nylon Sealed Ring

Not Included in SET

		
TWC-1105 Transfer For Wise Click	WCT-900 Wise Click Tool	AWC-1301 Analog For Wise Click




MULTI-UNIT STRAIGHT AESTHETIC ABUTMENTS


Abutments with smooth, non-locking surface for Multi-Unit Restoration

AESTHETIC REGULAR UNIT

Straight Multi-Unit

4.1 Regular Platform

		
ERU-4110 H = 1.5 mm	ERU-4120 H = 2.5 mm	ERU-4130 H = 3.5 mm



See the instructions Manual

MULTI SYSTEM

NOTE: The screw geometry was carefully designed and optimized to allow safe tightening and long service. Special attention must be paid while tightening the screw. Make sure to use an original SURCAM hex key with edges that do not show wear. It is recommended to purchase additional screws

MULTI-UNIT - ACCESSORIES

							
ISU-1201 H01	ISU-1202 H02	ISU-1203 H03	ISU-1801 H01	ISU-1802 H02	ISU-1803 H03	IBB-2530 Ball attachment	IBM-3230 Connection bar to be modified to the space between the sleeves
Ti Sleeve with 3 positioning options			Ti Sleeve with 2 positioning options				

							
PSU-0032 Plastic Sleeve screw included	TSU-1250 Titanium Sleeve screw included L = 12.0 mm	CMU-100006 Short slim titanium sleeve	CMU-100007 Long slim titanium sleeve L = 8.0 mm	PSU-1050 3D Scanning Lab device	PSU-100015 3D Scanning Intraoral device	PSU-0050 Parallel connection	TSM-100013 Titanium sleeve for welding

						
TMU-0532 Titanium Short Sleeve screw included	TMU-1000 Open Tray Short Transfer	TMU-1200 Open Tray Transfer screw included	TMU-8005 Closed-Tray Transfer	WQP-011 Plastic for Closed-Tray Transfer	HCC-4550 Healing Cap for Multi-Unit	AMU-1201 Analog for Multi-Unit

BALL ATTACHMENT ABUTMENTS

Overdenture Ball Attachment connection to Implants

BALL REGULAR STRAIGHT

4.1 Regular Platform Dia. 4.3 mm

				
BRSE-4120 H = 2.0 mm Without Cap	BRSE-4130 H = 3.0 mm Without Cap	BRSE-4140 H = 4.0 mm Without Cap	BRSE-4150 H = 5.0 mm Without Cap	BRSE-4160 H = 6.0 mm Without Cap

BALL REGULAR SET - NYLON

4.1 Regular Platform Dia. 4.3 mm

FULL SET

				
BRSE-4120-S Without Stainless Steel Housing & Nylon Insert	BRSE-4130-S Without Stainless Steel Housing & Nylon Insert	BRSE-4140-S Without Stainless Steel Housing & Nylon Insert	BRSE-4150-S Without Stainless Steel Housing & Nylon Insert	BRSE-4160-S Without Stainless Steel Housing & Nylon Insert

BALL HOUSING MINOR \ BALL NYLON INSERT

		
BHM-5032 Stainless Steel Housing	BHM-4730 Slim Stainless Steel Housing	BIS-0025 Nylon Insert

ANALOG FOR BALL ATTACHMENT


AMB-1201 Analog for Ball attachment

STRAIGHT AESTHETIC ABUTMENTS

Includes a Retaining Screw

In order to preserve the original structure of the titanium alloy, **SURCAM** utilizes advanced CAD/CAM methods in the production process. These methods effectively reduce undesirable deformations of the alloy structure.

SPECIFICATIONS

- The angular aesthetic Abutment body functions as impression coping
- Includes the Prosthetic Screw
- Installation using large external hex drivers (1.3mm)
- Non-rotational design for standard, internal-hex implants
- Compatible with ZIMMER®, MIS, ALPHA BIO, A.B.DENTAL, Bio Horizon
- Made of Commercial grade 5 Titanium alloy (ASTM-F-136 02A COND A)

STRAIGHT LEVEL EXTERNAL ABUTMENTS

4.1 Regular Platform

		
SLE-R09-10 H = 1.0 mm Straight Shoulder	SLE-R10-20 H = 2.0 mm Straight Shoulder	SLE-R11-30 H = 3.0 mm Straight Shoulder

15° ANGULAR CURVE EXTERNAL ABUTMENTS

4.1 Regular Platform

		
ACE-R15-10 H = 1.0 mm Aesthetic Shoulder	ACE-R15-20 H = 2.0 mm Aesthetic Shoulder	ACE-R15-30 H = 3.0 mm Aesthetic Shoulder

The Titanium / Plastic Abutment are for investment and casting following conventional crown and bridge techniques. The Titanium / Plastic Abutments combine a precision of a machined interface between the implant and the Abutment with the convenience of a Plastic sleeve for casting.

TITANIUM PLASTIC CASTING ABUTMENTS

4.1 Regular Platform

	
TPAE-1220 Hex Base Titanium Body	TPAE-1220-R Non-Engaging Titanium Body

5.0 Wide Platform

	
TPAWE-1220 Hex Base Titanium Body	TPAWE-1220-R Non-Engaging Titanium Body

Analogs are used for the fabrication of a master model.
They are replicas of the implant. All analogs are made of high grade stainless steel.

ANALOG FOR EXTERNAL HEX

4.1 Platform / 5.0 Platform

	
AES-1241 4.1 Platform Stainless Steel Body	AEWS-1251 5.0 Platform Stainless Steel Body

ANALOG FOR MULTI-UNIT


AMU-1201 Analog for Multi-Unit

ANALOG FOR WISE CLICK SYSTEM


AWC-1301 Analog for Wise Click

ANALOG FOR BALL ATTACHMENT


AMB-1201 Analog for Ball attachment

Healing caps are used to maintain access through tissue during the healing phase. They are available in a variety of heights and diameters to accommodate specific tissue depth.

HEALING CAP REGULAR PLATFORM

4.1 Platform

	
HCWE-5030 Dia. = 5.0 mm H = 3.0 mm	HCWE-5050 Dia. = 5.0 mm H = 5.0 mm

HEALING CAP WIDE PLATFORM


5.0 Platform

	
HCWE-6030 Dia. = 6.0 mm H = 3.0 mm	HCWE-6050 Dia. = 6.0 mm H = 5.0 mm

TRANSFERS IMPRESSION COPING

Includes a Retaining Screw

OPEN TRAY TRANSFER

	
TREQ-1101 Stainless Steel body 4.1 Platform	TWEO-1101 Stainless Steel body 5.0 Platform

TRANSFER WIDE CLOSED

Closed-tray Transfer Abutment

4.1 Platform

	
TCER-0910 Titanium Body 1.0 mm Aesthetic Shoulder	TCER-0930 Titanium Body 3.0 mm Aesthetic Shoulder

Plastic Cap Included

TRANSFER'S MATCHED PLASTIC

Closed-tray Transfer Plastic Cap

		
TMP-0008 Low Level Plastic Cap	TMP-0010 High Level Plastic Cap	CPC-0008 Burn Out Plastic Cap



This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.






Northern Industrial Zone P.O.B 12084
Nahariya 2201202, Israel
Technical support: 972-4-9523334
Sales: 972-4-9523511

www.surcamdental.com

E-mail: sales@surcamdental.com

 **Surcam Dental LTD**